

Orthodontic Burs

For de-bonding, adhesive removal and interproximal spacing

8 Flutes



US NO.	ISO NO.	1171	1172	169L	699L
FG (19MM)	500 314 194 006	012	016	009	009
RA (22MM)	500 204 194 006	012	016	009	
HP (44.5MM)	500 104 171 006			009	

Interproximal Diamond



FIG NO:	392
DIAMETER 1/10MM	016
FINE	•
MEDIUM	•

12 Flutes



US NO.	7006	7205	7642	7664	7901	7902	7903	OF1	274	48L	48L
FIG NO.	41	212	375R	375R	246	246	246	379	274	48L	48L
DIAMETER 1/10MM	018	016	010	014	009	010	012	023	016	010	012
H/LENGTH (MM)	1.8	8.5	7.0	8.2	3.8	3.8	3.9	4.3	3.5	8.0	8.0

12 Flutes Safe End





FIG NO:	H379-023AGK	H22-016ALGK
FG	023	016
RA	023	016

18 Flutes



DIAMETER 1/10MM	016 (SHORT)	016 (REG.)	016 (LONG)
INTERNATIONAL REF.	118S	218	118L
HEAD LENGTH (MM)	5.7	7.3	8.3
OVERALL LENGTH (MM)	19	22	22